# Withings BPM Connect

## INSTRUCTION MANUAL

---

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EN | FR | DE

# THANK YOU FOR CHOOSING WITHINGS BPM CONNECT

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This guide is available in additional languages at: [www.withings.com/guides](http://www.withings.com/guides)

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
## 1 Visit [go.withings.com](https://go.withings.com) from your smartphone or tablet and follow the instructions to install the Health Mate app.

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## 2 In the app, follow the instructions to install Withings BPM Connect.

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# How to take a blood pressure measurement

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
2

**Unroll cuff. Place your left arm inside it.**

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## 3

**Tighten the cuff around your arm. The tube should be positioned against your inner arm.**

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
4

**Place your arm on a table and level with your heart.**

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## 5

**Turn on your blood pressure monitor to start the measurement. You can select a single or triple measurement with a long press on the button. Press the button again to start the measurement.**

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
10

## 6

**If ERROR message appears on screen, please refer to the explanation displayed in the app and troubleshooting described from page 15 onward.**

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
## Instructions manual

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## Requirements

To install your BPM Connect, you will need a device with iOS (10.0 or higher) or Android (6.0 and higher). After installation, you do not need your phone to take a measurement. The device will sync via Wi-Fi or Bluetooth\*.

## Important safety information

Consult your doctor during pregnancy, or if you suffer from arrhythmia and arteriosclerosis. Please read this section carefully before using the Wi-Fi Smart Blood Pressure monitor.

## Intended use

The BPM Connect is a digital monitor intended for use in measuring blood pressure and heart rate. The device is intended to be used in a human adult population with an arm circumference between 9 inches to 17 inches (22 cm to 42 cm).

The BPM Connect is a medical device. The BPM Connect is intended to be used at home. The BPM Connect must be activated with the Withings Health Mate companion app. Withings Health Mate companion app is an accessory to the Withings BPM Connect and is not a medical device. If you observe hypertensive values, contact your physician.

## General safety and precautions

- Do not forcibly bend the arm cuff.
- Do not inflate the arm cuff when it is not wrapped around your arm.
- Do not drop or apply strong shocks and vibrations to the blood pressure monitor.
- Do not take measurements after bathing, drinking alcohol, smoking, exercising or eating.
- Do not immerse the arm cuff in water.
- Do not use with a pacemaker, defibrillator or other electric implant.
- Use the device on adults only.
- Do not use on children or pets.

Android is a trademark of Google Inc.  
iOS is a trademark of Apple Inc.

## General usage

Always consult your doctor if you believe you have a blood pressure disorder. People with severe blood flow problems, or blood disorders, should consult a doctor before using the blood pressure monitor. Cuff inflation may cause internal bleeding. Conditions such as common arrhythmias, ventricular premature beats, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, pre-eclampsia or renal disease may affect the performance of the automated sphygmomanometer and/or its blood pressure reading. BPM Connect is a precision measuring device designed to be understood by lay users but it should still be handled with care. Exposing the device to lint, dust or sunlight for a long period of time may reduce its lifetime or damage it. A damaged cuff or sensor may lead to incorrect measurements. The device may not be connected via a USB port. This device is intended for home use and is not to be used in a clinical setting. The cuff will come in contact with the skin. Measurements may be affected by extreme temperatures, humidity and altitude.

Following REACH Regulation 1907/2006, there is an open risk on plumb and D6.

- Do not leave the blood pressure monitor unattended with infants or persons who cannot express their consent.
- Do not use the blood pressure monitor for any purpose other than measuring blood pressure.
- Do not disassemble the blood pressure monitor.
- Do not operate the blood pressure monitor in a moving vehicle such as a car or an airplane.
- Do not use a cellular phone near the device.
- Do not use the device with the USB cable plugged.
- Using the cuff or too frequent measurements continuously or taking measurements too frequently may interfere with blood flow and result in injury to the user.
- Do not apply the cuff if the arm has a wound or is undergoing medical treatment, as this can cause further injury.
- The cuff should not be used on the arm on the side of a mastectomy.
- Using of the device with other equipment on the same limb may cause damage to the device or other equipment.

## **Cleaning**

- Do not use an alcohol-based cleaning solution or solvent to clean the device
- Clean the device with a soft, dry cloth
- If the cuff is dirty, clean it with a moistened cloth and mild soap.
- Do not flush the device and cuff with water.
- Do not dismantle the device, disconnect the cuff, or try to repair it yourself. If the device is not working, contact Withings support.
- Do not operate the device in extreme temperatures, humidity, or direct sunshine
- Do not violently shake the unit or drop it on the floor.
- Do not submerge the device or any of the components in water

## **Storage**

- Store the device and the components in a clean and safe location
- If storage conditions are different from the usage conditions indicated in this document, please wait 30 minutes before taking a measurement.

## **Maintenance**

If you cannot fix a problem using the troubleshooting instructions, please request service from your dealer or Withings support. The manufacturer will make available on requested circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist staff or authorized representatives with repair.

It is generally recommended to have the device inspected every 2 years, to ensure proper functioning and accuracy. Do not use the device while it is under maintenance.

---

Error message	Causes	Countermeasures
<b>The measurement could not be performed. Please try again.</b>	Hardware auto-check error.	The measurement could not be performed. Please try again. If the problem occurs again, please contact customer service.
<b>Connection error.</b>	There is a communication error between the blood pressure monitor and the collateral device.	Push the blood pressure monitor button for 3 seconds to turn it off. Quit the application, turn your bluetooth connection off and on. Push the blood pressure monitor button to turn it on.
<b>Blood pressure measurement could not be performed</b>	<ul style="list-style-type: none"><li>- Cuff is not deflated before measurement.</li><li>- Noise interference.</li><li>- Excessive user movements.</li></ul>	Please wait until the cuff is totally deflated before starting a new measurement. Stay still and do not speak during the measurement.

---

---

Error message	Causes	Countermeasures
<b>Check that the blood pressure monitor is correctly positioned on your arm and that the measurement is being performed under the proper conditions.</b>	The measurement may not occur because of: <ul style="list-style-type: none"><li>- Broken cuff</li><li>- Pump or valve failure</li><li>- Pressure overflow</li><li>- Noise interference</li><li>- User movement during the measurement</li><li>- Special characteristics, such as existing heart issues</li><li>- Low battery level</li></ul>	<ol style="list-style-type: none"><li><b>1.</b> Check that the Blood Pressure Monitor is correctly positioned on your arm and that measurement is performed in good conditions.</li><li><b>2.</b> Low battery level. Load your BPM Connect with a micro USB cable.</li><li><b>3.</b> If user has special characteristics, please contact your physician.</li><li><b>4.</b> If the problem continues to occur, please contact customer service.</li></ol>
<b>Low battery level.</b>	Battery level is too low.	Load the device with the provided charging cable.

---

---

Problem	Remedy
<b>The cuff cannot be inflated or the air pressure cannot rise</b>	<ol style="list-style-type: none"><li data-bbox="720 273 1684 296">1. Check the cuff position, fasten the cuff correctly and remeasure the pressure again</li><li data-bbox="720 304 1684 327">2. Check the blood pressure monitor connection to the iOS or Android device</li></ol>
<b>The low battery icon is displayed on the device</b>	Reload the device with the cable provided
<b>The measurement cannot be taken and the application shows an error message</b>	<ol style="list-style-type: none"><li data-bbox="720 449 934 472">1. Re-fasten the cuff</li><li data-bbox="720 480 953 503">2. Relax and sit down</li><li data-bbox="720 511 1478 534">3. Keep the cuff level with your heart during the measurement period</li><li data-bbox="720 542 1286 565">4. Keep silent and stay still during the measurement</li><li data-bbox="720 573 1633 596">5. If you have severe heart issue, the blood pressure may not be measured correctly</li></ol>
<b>Under normal circumstances, the measurement may differ from the measurement taken in a medical setting. Your blood pressure may change because of the environment in which the measurement is taken. Please share the values recorded via BPM Connect with your healthcare practitioner.</b>	<ol style="list-style-type: none"><li data-bbox="720 646 1305 669">1. The variation is due to the different environments</li><li data-bbox="720 677 1755 731">2. The blood pressure is changing according to the physiological or psychological status of the human body</li><li data-bbox="720 739 1258 762">3. Show your recorded values to your physician</li></ol>

---

# Specifications Technical data

<b>Product description:</b>	Digital automatic Blood Pressure Monitor
<b>Model:</b>	WPM05
<b>Blood pressure measurement method:</b>	Cuff oscillometric method
<b>Cuff inflation:</b>	Automatic inflation with air pump at 15 mmHg/s
<b>Pressure sensor:</b>	Gauge sensor
<b>Measurement range (pressure):</b>	0 to 285 mmHg, DIA 40 to 130 mmHg, SYS 60 to 230 mmHg
<b>Measurement range (pulse):</b>	40 to 180 beats/min
<b>Accuracy (pressure):</b>	Within +- 3 mmHg or 2% of reading
<b>Accuracy (pulse):</b>	Within +-5% of reading
<b>Sensor:</b>	Semiconductor pressure sensor
<b>Operating conditions:</b>	5 to 40°C, 15 to 90% RH, atmospheric 86Kpa-106kpa, altitude: 2000m
<b>Storage and transport conditions:</b>	-20 to 60°C, 15 to 93% RH, atmospheric 50Kpa-106kpa, altitude: 2000m

<b>Arm type:</b>	Use on left arm.
<b>Power source:</b>	5V 1A
<b>Weight:</b>	Approx. 245g
<b>Accessories:</b>	Charging cable, Instruction manual
<b>Product life:</b>	3 years
<b>Wireless transmission:</b>	Wi-Fi and BLE

**Note:** Specifications are subject to change without prior notice or any obligation on the part of the manufacturer

**Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal.**

**The symbol applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste.**

At the end of the device's useful life, the user must deliver it to a collection center for electric and electronic garbage, or return it to the retailer when purchasing a new device. Disposing of the product separately prevents possible negative environmental and health consequences deriving from inadequate disposal. It also allows the recovery of materials it is composed of to save energy and resources and avoid negative effects to the environment and health. In case of abusive disposal of device by the user, will be applied administrative endorsements in compliance with current standards. The device and its parts must be disposed of as appropriate, in accordance with national or regional regulations.

Version 1.0 April 2019

## Federal Communications Commission (FCC) Statement 15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

### 15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance.

**This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. This device complies with Industry Canada license-exempt RSS standards. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.**

**FCC ID: XNAWPM05**

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following.

| Interference may occur in the vicinity of equipment marked with | Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. | The use of accessories and cables other than those specified may result in increased emissions or decreased immunity | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter | The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Any other accessories, transducers and cables may result in increased emissions or decreased immunity and EMC performance.

| The device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, which should be observed to verify normal operation in the configuration in which it will be used. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following. Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

## Guidance and manufacturer's declaration-electromagnetic emissions

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The Wi-Fi Smart Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RE emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The Wi-Fi Smart Blood Pressure Monitor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

# Declaration – electromagnetic emissions and immunity 23

for equipment and systems that are not life-supporting and are specified for use only in a shielded location

---

## The Wi-Fi Smart Blood Pressure Monitor declaration electromagnetic immunity

---

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should ensure that it is used in such an environment.

---

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RFIEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A N/A	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT or SYSTEM including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol.

---

# Declaration – electromagnetic emissions and immunity

for equipment and systems that are not life-supporting and are specified for use only in a shielded location

---

## Declaration – electromagnetic immunity

---

The Wi-Fi Smart Blood Pressure Monitor is intended for use in the electromagnetic environment specified below.

The customer or the user of the Wi-Fi Smart Blood Pressure Monitor should assure that it is used in such an environment.

---

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: $\pm 8$ kV Air: $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	Contact: $\pm 8$ kV Air: $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): $\pm 1$ kV line(s) to earth: $\pm 2$ kV 100 kHz repetition frequency	line(s) to line(s): $\pm 1$ kV 100 kHz repetition frequency	A Mains power quality should be that of a typical commercial or hospital environment.

---

---

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT for 0.5 cycle, At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT ; 25/30 cycles Single phase at 0° 0% UT; 300 cycle	0% UT for 0.5 cycle, At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT ; 25/30 cycles Single phase at 0° 0% UT; 300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EQUIPMENT or SYSTEM requires continued operation during power mains interruptions, it is recommended that the EQUIPMENT or SYSTEM be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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## **Withings Two (2) Years commercial guarantee Withings BPM Connect | Wi-Fi Smart Blood Pressure Monitor**

Withings SA, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, France ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

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**This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s).**

**Operation is subject to the following two conditions:**

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

**Radiation Exposure Statement:** The product comply with the Canada portable RF exposure limit set forth for an uncontrolled environment and are safe for intended operation as described in this manual. The further RF exposure reduction can be achieved if the product can be kept as far as possible from the user body or set the device to lower output power if such function is available.

This equipment should be installed and operated with minimum distance 0cm between the radiator & your body.

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**This device complies with the following normative documents : COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 concerning medical devices as amended by Directive 2007/47/EC**

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EN ISO 13485: Medical devices - Quality management systems - Requirements for regulatory purposes. Reference to standards contd. | EN ISO 14971: Medical devices - Application of risk management to medical devices. | IEC/EN 60601-1: Medical electrical equipment - Part 1: General requirements for basic safety and essential performance. | EN 1060-3: Non-invasive sphygmomanometers, Part 3: Supplementary requirements for electromechanical blood pressure measuring systems. | EN 1060-4: Non-invasive sphygmomanometers. Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers. | IEC/EN 60601-1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. | IEC 80601-2-30: Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers. | EN 300 328: Electromagnetic compatibility

and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the Directive (2014/53/EU). | EN 301 489-1: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements. | EN 301 489-17: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements | IEC/EN60601-1-2: Medical electrical equipment: Part 1-2: General requirements for basic safety and essential performance-collateral standard electromagnetic compatibility | EN 55011: Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement. | FCC part B 15B: Electromagnetic Compatibility. | FCC Rule Part: 15.247 Cat: DSS (Bluetooth) FCC Rule Part: 15.247 Cat: DTS (BT4.0). | EN ISO 10993-1: Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process. | EN ISO 10993-5: Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity | EN ISO 10993-10: Biological evaluation of medical devices. Tests for irritation and skin sensitization | IEC/EN 60601-1-6: Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability | ANSI/AAMI/ISO 81060-2 Non-invasive sphygmomanometers Part 2: Clinical validation of automated measurement type | EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

**Distributed by**

Withings Inc.,  
179 South Street, Floor 5  
Boston, MA 02111 USA

2019

**Manufacturer/Fabricant/Hersteller:**

WITHINGS SA,  
2 rue Maurice Hartmann,  
92130 Issy-les-Moulineaux, FRANCE  
+33141460460

**IP22****Ingress of water or particulate matter**

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**Complies with waste electrical and electronic equipment directive**

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Frequency band :  
2402 - 2480 MHz  
Maximum output power : 4.45 dBm

**EN** Withings SA hereby declares that the device Withings BPM Connect is in conformity with the essential requirements and other relevant requirements of Directive 2014/53/EU.

The declaration of conformity can be found at: [withings.com/compliance](http://withings.com/compliance).

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Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Australian sponsor:

Emwergo Australia

Level 20, Tower II, Darling Park, 201 Sussex Street, Sydney, NSW 2000  
Australia

## Need help?

Besoin d'aide? | Brauchen Sie Hilfe?

➔ [withings.com/support](https://withings.com/support)

withings

Withings BPM Connect | Wi-Fi Smart Blood Pressure Monitor | Apr.2019

© 2019 Withings SA. All rights reserved.