
Withings Sleep Analyzer

QUICK INSTALLATION GUIDE

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EN | LV | LT | ET

THANK YOU FOR CHOOSING WITHINGS SLEEP ANALYZER

PALDIES, KA IZVĒLĒJĀTIES WITHINGS SLEEP ANALYZER | DÉKOJAME, KAD PASIRINKOTE WITHINGS SLEEP ANALYZER |
AITĀH, ET VALISITE SEADME WITHINGS SLEEP ANALYSER.

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

USB Power Adapter

LV USB strāvas adapteris
LT USB maitinimo adapteris
ET USB-toiteadapter


How to install your Withings Sleep Analyzer

Kā uzstādīt Withings Sleep Analyzer | Kaip pradėti naudoti Withings Sleep Analyzer |
Uneanalüsaatori Withings Sleep Analyzer paigaldamine

1

Place the entire sleep sensor under your mattress, so that it will be positioned underneath your chest. Only the cable should extend beyond the edge of your mattress.

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
2

Plug the sensor into an outlet using the provided adapter.

- LV Pievienojiet sensoru kontaktligzdai, izmantojot komplektācijā iekļauto adapteri.
- LT Prijunkite jutiklį prie lizdo naudodami pateiktą adapterį.
- ET Ühendage andur komplekti kuuluva adapteri abil pistikupessa.

How to install your Withings Sleep Analyzer

Kā uzstādīt Withings Sleep Analyzer | Kaip pradėti naudoti Withings Sleep Analyzer |
Uneanalüsaatori Withings Sleep Analyzer paigaldamine

3

Visit go.withings.com from your smartphone or tablet, and follow the instructions to install the Health Mate app.

- LV** Apmeklējiet tīmekļa vietni go.withings.com viedtālrunī vai planšetdatorā un izpildiet norādījumus, lai instalētu lietotni Health Mate.
- LT** Savo išmaniajame telefone arba planšetiniame kompiuteryje apsilankykite go.withings.com svetainėje ir sekite „Health Mate“ programėlės įdiegimo instrukcijas.
- ET** Külstage nutitefonis või tahvelarvutis veebisaiti go.withings.com ja järgige rakenduse Health Mate installimisjuhiseid.

4

In the Health Mate app, follow the instructions to install Withings Sleep Analyzer.

- LV Lietotnē Health Mate izpildiet norādījumus, lai instalētu Withings Sleep Analyzer.
- LT Health Mate programēlėje sekite instrukcijas, kad įdiegtumėte Withings Sleep Analyzer.
- ET Avage rakendus Health Mate ja järgige Withings Sleep Analyzeri installimisjuhiseid.

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
To clean the textile cover of the sensor, gently remove the cover first. The cover can be machine washed, but do not tumble dry or iron it. Wait until the cover is completely dry before replacing it on the sensor.

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
LT Sąrankos LED

ET Seadistus-LED

EN You can deactivate the sound sensor in your settings in the Health Mate app

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Federal Communication Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Important note: Radiation Exposure Statement

The product complies with the US portable radio frequency (RF) exposure limit set forth for an uncontrolled environment and is safe for intended operation as described in this manual. Further RF exposure reduction can be achieved if the product is kept as far as possible from the user's body or if the device is set to lower output power if this function is available.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Country Code selection feature to be disabled for products marketed to the US.

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IC: 11411A-WSM02

This device complies with Industry Canada's licence exempt RSSs. Operation is subject to the following two conditions:

- (1) This device may not cause interference; and
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

Radiation Exposure Statement

The product complies with the Canada portable RF exposure limit set forth for an uncontrolled environment and is safe for intended operation as described in this manual. Further RF exposure reduction can be achieved if the product is kept as far as possible from the user's body or if the device is set to lower output power if such function is available.

Country Code selection feature to be disabled for products marketed to CANADA.

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US Withings Two (2) Years Limited Warranty - Withings Sleep Analyzer

Withings SA 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, France ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

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UK Withings Two (2) Years commercial guarantee - Withings Sleep Analyzer

Withings SA 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, France ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

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© 2020 Withings SA. All rights reserved.