
Withings BPM Core

INSTRUCTION MANUAL

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EN | FR | DE

THANK YOU FOR CHOOSING WITHINGS BPM CORE

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This guide is available in additional languages at: www.withings.com/guides

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
1

Visit go.withings.com from your smartphone or tablet and follow the instructions to install the Health Mate app.

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
2

In the app, follow the instructions to install Withings BPM Core.

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8-9

1

BPM Core allows you to take two types of measurements:

- **Blood pressure, which uses the cuff as a sensor**
- **ECG and stethoscope signal recording, which uses three electrodes (two inside the cuff, the other on the metal tube) and the stethoscope (the circular white sensor on the cuff)**

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Before the measurement

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1

Before the measurement, make sure that you:

- Use the product on the left upper arm.
- Rest for 5 minutes before the measurement.
- Sit down in a comfortable position, legs uncrossed. feet flat on the floor, arm and back supported.
- Do not speak or move during the measurement.
- Wear only one layer of clothing, which should not cover your left arm. The electrodes should be in contact with the skin.
- Take the measurement in a calm and quiet area.

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
12-13

1

Tighten the cuff around your arm. The stethoscope should be facing the rib cage, with the metal tube facing down about one inch above the elbow.

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Left arm position depends on the type of measurement being performed

1

For blood pressure measurement

Place your left arm on a table with the cuff level with your heart and your palm facing up. Make sure the cuff is not touching the left side of your chest. Your arm should be resting comfortably on the table without making any effort. Do not speak or move during the measurement.

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
Left arm position depends on the type of measurement being performed

2

For ECG and stethoscope signal recording :

The excess cuff should not overlap the stethoscope. For this measurement, you should place the stethoscope in contact with your chest, and the two electrodes inside the cuff should be in contact with your skin. Make sure that the stethoscope is in contact with your chest during the measurement. You should not move or talk during the measurement. Breathe normally. You do not need to put your left hand on a table.

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
Mode standard | Standard-Einstellung

1

Two measurements in a row will be taken during this mode. The first one is a blood pressure measurement, which you should take in the position described on p.17 The second is a combined measurement of an ECG and heart sounds check via the stethoscope.

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2

Push the button. The screen will display "START." Push the button once again to start the measurement.

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Mode standard | Standard-Einstellung

3


When the screen displays the illustrated pictogram, it means that the blood pressure measurement is finished and that you will now proceed to the ECG and stethoscope signal recording.

Adjust your arm position to place the stethoscope against your chest, and place your hand on the metal tube to ensure an accurate measurement. You can remove your hand from the table and follow the position described on the p.17 .

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1

When you push the button to start BPM Core, you can also select “Start x3” or “ECG” by swiping the light trail on the screen. “Start x3” launches three blood pressure measurements in a row followed by an ECG and stethoscope measurement. “ECG” will only record the ECG and stethoscope signals. Make sure to use the correct arm position for each measurement type.

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The LED of the button of BPM Core can be:

Blue: if BPM Core has not been paired with a mobile device yet.

Green: for optimal and normal blood pressure measurements. The LED can also be green if BPM Core has been paired with a mobile device.

Orange: for high normal blood pressure measurements. The LED can also be orange if BPM Core has a low battery (you can still take a measurement) or if there is an error during one of the measurements during x3 measurement mode. During x3 mode, you may see a short ERROR display and orange LED, but afterward, it will give you the option to restart the measurement. If you fail to get a result on the second try, the LED will turn red.

Red: if hypertension (grade 1, 2, 3) or isolated systolic hypertension are detected. The LED can also be red if BPM Core has no battery or if there is an error during the blood pressure or stethoscope measurement.

Pink: if you press and hold the monitor's button for 6 seconds until the Settings menu opens and "Training" is displayed on the screen. The LED can also be pink when BPM Core has been factory reset.

Note: The LED color displayed on the device may vary depending on the classification guidelines of your country. Please go to support.withings.com for more information on blood pressure measurement classification.

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
Message d'erreur | Fehlermeldung

26-27

If ERROR message appears on screen, please refer to the explanation displayed in the app and troubleshooting procedures described beginning on p. 34.

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
30-31

Requirements

In order to use your BPM Core, you need an iOS (10.0 or higher) or Android (6.0 and higher) device to install it. Thereafter, the product can be used without your mobile device on you, thanks to the Wi-Fi and Bluetooth® connections.

Important safety information

Consult your doctor during pregnancy, or if you suffer from arrhythmia and arteriosclerosis. Please read this section carefully before using this smart blood pressure monitor with ECG & digital stethoscope.

Intended use

BPM Core is a digital monitor intended for use in measuring blood pressure and heart rate, to record a one-lead ECG to detect atrial fibrillation, and to uncover valvular heart diseases (aortic stenosis, aortic regurgitation, mitral regurgitation). The device is intended to be used in a

human adult population with an arm circumference of 9 inches to 17 inches (22 cm to 42 cm). BPM Core is a medical device. Contact your physician if hypertensive values, atrial fibrillation (AFib) or VHD are indicated.

General safety and precautions

- Do not forcibly bend the arm cuff.
- Do not inflate the arm cuff when it is not wrapped around your arm.
- Do not apply strong shocks and vibrations to the blood pressure monitor or drop it.
- Do not take measurement after bathing, drinking alcohol, smoking, exercising or eating.
- Do not immerse the arm cuff in water.
- Do not use with a pacemaker, a defibrillator or other electric implant.
- Use on adults only.
- Do not use on children or pets.

Android is a trademark of Google Inc. iOS is a trademark of Apple Inc.

Caution

EN

General usage

Always consult your doctor. Self-diagnosis of measurement results and self treatment are dangerous. People with severe blood flow problems or blood disorders should consult a doctor before using the blood pressure monitor. Cuff inflation can cause internal bleeding. Operational factors such as common arrhythmias, ventricular premature beats, arterial sclerosis, poor perfusion, diabetes, age, pregnancy, pre-eclampsia or renal disease can affect the performance of the automated sphygmomanometer and/or its blood pressure reading. This is a precision measuring device that may be understood by lay users, but should still be handled with care. Exposing the device to prolonged lint, dust or sunlight might reduce its life or damage it. A damaged cuff or sensor may lead to incorrect measurements. The

USB port should only be used for charging the device. The user is an intended operator. Parts in contact with the skin: cuff and electrodes. Measurements can be affected by extreme temperatures, humidity & altitude.

Caution

EN

32-33

- Do not leave the blood pressure monitor unattended with infants or people who cannot express their consent.
- Do not use the blood pressure monitor for any purpose other than measuring blood pressure, recording an ECG and recording heart sounds.
- Do not disassemble the blood pressure monitor.
- Do not operate the blood pressure monitor in a moving vehicle (car, airplane).
- Do not use the device with the USB cable plugged in.
- Improper continuous pressure of cuff or too-frequent measurements may interfere with blood flow and result in injury to the user. Check to ensure that the use of the device does not result in prolonged impairment of your blood circulation.
- Do not apply the cuff over the user's arm if it has a wound or medical treatment, as this can cause further injury.
- Consult your physician before using this monitor if you have had a mastectomy.
- Use of the device can temporarily cause the loss of

function of equipment that is used simultaneously on the same limb.

- Do not use if you have a congenital cardiomyopathy, aorto-venous fistula, pulmonary pathology, a prosthetic heart valve, if you are on kidney dialysis, or if you are hemodynamically unstable.
- There is a risk concerning Lead, PZT, DEHP, DBP, SCCP and D6 according to REACH regulation 1907/2006.

After use

EN

Cleaning

- Do not use an alcohol-based or solvent agent to clean the device.
- Clean the device with a soft and dry cloth.
- The dirt on the cuff can be cleaned with a damp cloth and soap.
- Do not flush the device and cuff with a lot of water.
- Do not dismantle the device, disconnect the cuff, or try to repair it by yourself. If a problem occurs, please check with the distributor.
- Do not operate the device in a severe environment of extreme temperature, humidity, or direct sunshine.
- Do not shake the unit violently.
- Do not submerge the device or any of the components in water.
- Do not use the device after a strong shock, such as dropping the unit on the floor.

Storage

- Store the device and the components in a clean and safe location.
- If storage conditions are different from the usage conditions indicated in this document, please wait 30 minutes before taking a measurement.

Maintenance

If you cannot fix the problem using the troubleshooting instructions, request service from your dealer. The manufacturer will make available on request circuit diagrams, component part lists, descriptions, calibration instructions, or other information that will assist the manufacturer's staff or authorized representative with repair. It is generally recommended to have the device inspected every 2 years to ensure proper functioning and accuracy. Do not use the device while doing maintenance steps.

Error messages

EN

34-35

Error message	Causes	Solutions
Measurement could not be performed. Please try again. If the problem occurs again, please contact customer service.	Hardware auto-check error	Measurement could not be performed. Please try again. If the problem occurs again, please contact customer service.
Connection error. Quit the application and reconnect your blood pressure monitor. Connection error.	Communication error between the blood pressure monitor and the collateral device.	<ol style="list-style-type: none">1. Push the blood pressure monitor button for 3 seconds to turn it off. Quit the application and push the blood pressure monitor button to turn it on.2. If the problem occurs again, please contact customer service.
Please wait until the cuff is totally deflated before starting a new measurement. Stay still during measurement.	<ul style="list-style-type: none">- Cuff is not deflated before measurement- Noise interference- Excessive user movements	<ol style="list-style-type: none">1. Please wait until the cuff is totally deflated before starting a new measurement. Stay still during the measurement2. If the problem occurs again, please contact customer service.

Error messages

EN

Error message	Causes	Solutions
Check that the blood pressure monitor is correctly positioned on your arm and that the measurement is being performed under good conditions.	<ul style="list-style-type: none">- User movement, or user is not relaxed, or is talking- Cuff is not correctly fastened- Cuff is broken- Pump or valve failure- Pressure overflow- Noise interference- User has special characteristics, e.g., some forms of arrhythmia can trigger errors while measuring blood pressure- Battery level is too low	<ol style="list-style-type: none">1. Check that the blood pressure monitor is correctly positioned on your arm and that the measurement is being performed under good conditions.2. Low battery level. Charge your BPM Core with a micro USB cable.3. If user has special characteristics, please contact your physician.4. If the problem occurs again, please contact customer service.
Low battery level	Battery level is too low	Charge the device with the provided charging cable

Troubleshooting

EN

36-37

Problem	Troubleshooting
No results for the ECG measurement	Ensure that the three electrodes are in direct contact with the skin. If yes, clean them with a tissue
A valvular heart disease analysis cannot be performed	Check the position of the stethoscope. The flat surface must be in contact with the chest. Ensure that you performed the required number of measurements to get a reading.
The inflation action cannot be performed or the air pressure cannot rise	<ol style="list-style-type: none">1. Check the cuff position, fasten the cuff correctly and remeasure blood pressure again2. Check the blood pressure monitor's connection to the iOS or Android device
The low battery icon is displayed on the device	Charge the device with the cable provided
The blood pressure cannot be taken and the application shows an error message	<ol style="list-style-type: none">1. Re-fasten the cuff2. Relax and sit down3. Keep the cuff at the same level as your heart during the measurement period4. Keep silent and still during the measurement5. If you have a severe heartbeat problem, the blood pressure may not be read correctly
Under normal circumstances, the measurement taken at home may differ from the measurement taken in a medical setting	<ol style="list-style-type: none">1. The variation is due to the different environments2. The blood pressure is changing according to the physiological or psychological status of the person being measured3. Atrial fibrillation may not be permanent in its early stages4. Show your recorded values to your physician

Specifications: technical data

EN

Product description	Automatic electronic blood pressure monitor with stethoscope, electrocardiogram and heart rate sensor
Model	WPM04
Blood pressure measurement method	Cuff oscillometric method
Cuff inflation	Automatic inflation with air pump at 5 mmHg/s
Pressure sensor	Gauge sensor
Measurement range (pressure)	0 to 285 mmHg, DIA 40 to 130 mmHg, SYS 60 to 230 mmHg
Measurement range (pulse)	40 to 180 beats/min
Pressure sensor accuracy	Within +- 3 mmHg or 2% of reading
Clinical accuracy (blood pressure)	BPM Core performance in blood pressure measurement was validated in a clinical trial. The results are within the margin of acceptance defined by the internationally recognized evaluation standard of blood pressure monitors ANSI/AAMI/ISO 81060-2:2013, EN ISO 81060-2:2014, developed by the European Society of Hypertension, British Hypertension Society and Association for the Advancement of Medical Instrumentation/American Heart Association. For more information, go to: withings.com/support

Specifications: technical data

EN

38-39

Accuracy (pulse)	Within +-5% of reading
Sensor	Semiconductor pressure sensor
Operating conditions:	5 to 40°C, 15 to 90% RH, atmospheric 86Kpa-106kpa, altitude: 2000m
Storage and transport conditions:	-25 to 70°C, 10 to 95% RH, atmospheric 86Kpa-106kpa, altitude: 2000m
Arm type	Use on left arm
Power source	3.6 VDC Lithium ion battery
Weight	Approx. 430g
Accessories	Charging cable, instruction manual
Product life	3 years
Wireless transmission	Wi-Fi and BLE

Note: Specifications are subject to change without prior notice or any obligation on the parts of the manufacturer

Disposal

EN

Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal. The symbol applied on the device or its packaging means that at the end of its useful life, the product must not be disposed of with domestic waste.

At the end of the device's useful life, the user must deliver it to collection centers for electric and electronic garbage, or give it back to the retailer when purchasing a new device. Disposing of the product separately prevents possible negative consequences for the environment and for health deriving from inadequate disposal. It also allows the recovery of materials of which the device is made up in order to obtain an important saving of energy and resources and to avoid negative effects to the environment and health. Failure to comply with the rules for sorting or recycling waste may make the user liable for a fine. The device and its parts are made with regard to disposal, as appropriate, in accordance with national or regional regulations.

Version 1.1 September 2020

FCC Statement

EN

40-41

Federal Communications Commission (FCC) Statement 15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.

- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1) this device may not cause harmful interference and 2) this device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. End users must follow the specific operating instructions for satisfying RF exposure compliance.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

RF Statement

EN

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following.

Interference may occur in the vicinity of equipment marked with the following symbol: | Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. | The use of accessories and cables other than those specified may result in increased emissions or decreased immunity | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes | Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than

the recommended separation distance calculated from the equation applicable to the frequency of the transmitter | The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Any other accessories, transducers and cables may result in increased emissions or decreased immunity and EMC performance. | The device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, which should be observed to verify normal operation in the configuration in which it will be used. Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following. Portable and mobile RF communication equipment (e.g. cell phones) can affect Medical Electrical Equipment. The use of accessories and cables other than those specified may result in increased emissions or decreased immunity of the unit.

RF Statement

EN

42-43

Guidance and manufacturer's declaration – electromagnetic emissions

This smart blood pressure monitor with ECG & digital stethoscope is intended for use in the electromagnetic environment specified below. The customer or the user of the smart blood pressure monitor with ECG & digital stethoscope should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
CE emissions CISPR11	Group 1	This smart blood pressure monitor with ECG & digital stethoscope uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RE emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	This smart blood pressure monitor with ECG & digital stethoscope is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	

Declaration — electromagnetic emissions and immunity

for equipment and systems that are not life-supporting and are specified for use only in a shielded locationshielded location

Smart blood pressure monitor with ECG & digital stethoscope declaration of electromagnetic immunity

This smart blood pressure monitor with ECG & digital stethoscope system is intended for use in the electromagnetic environment specified below. The customer or the user of the smart blood pressure monitor with ECG & digital stethoscope system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Conducted RFIEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	N/A
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the equipment or system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol. 

Declaration – electromagnetic immunity

This smart blood pressure monitor with ECG & digital stethoscope system is intended for use in the electromagnetic environment specified below. The customer or the user of the smart blood pressure monitor with ECG & digital stethoscope system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A N/A	The main power quality should be similar to that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	N/A N/A	The main power quality should be similar to that of a typical commercial or hospital environment.

Declaration — electromagnetic emissions and immunity

for equipment and systems that are not life-supporting and are specified for use only in a shielded location shielded location

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	-5% UT (95% dip in UT) for 0.5 cycle, -40% UT (60% dip in UT) for 5 cycles, -70% UT (30% dip in UT) for 25 cycles, -5% UT (95% dip in UT) for 5 sec	N/A	The main power quality should be similar to that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power main interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ECG performance

EN

46-47

Clinical performance

The ECG algorithm ability to accurately classify an ECG into AFib and sinus rhythm was tested in a clinical trial of 115 subjects. Rhythm classification of a 12-lead ECG by a cardiologist was compared to the simultaneous classification of an ECG recorded by BPM Core.

	Sensitivity(%)	Specificity
Atrial fibrillation	94.57	93.56

	Sensitivity (%)	Positive predictivity (%)
Normal sinus rhythm	93.76	93.78

The ECG recording result on the mobile app provides a detailed display of the result.

BPM Core outputs

Measurement	Display	Interpretation
Blood pressure	SYS 127	Systolic blood pressure value in mmHG
	DIA 82	Diastolic blood pressure value in mmHG
	HR 65	Heart rate value in beats per minute
ECG	ECG NORMAL	ECG signal shows a normal sinus rhythm
	ECG AFIB	ECG signal shows a risk of atrial fibrillation
	ECG INCONCLUSIVE	ECG signal is inconclusive and can't be interpreted by the algorithm

Information related to VHD detection is available in the Health Mate app.

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Garantie | Garantie

89

UK | Withings Two (2) Years Commercial Guarantee - Withings BPM Core | Smart Blood Pressure Monitor with ECG & Digital Stethoscope

Withings, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, France ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or Quick Installation Guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

US | Withings One (1) Year Limited Warranty - Withings BPM Core | Smart Blood Pressure Monitor with ECG & Digital Stethoscope

Withings, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of ONE (1) YEAR from the date of original retail purchase by the end-user purchaser ("Warranty Period"). Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or Quick Installation Guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use.

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This device complies with Part 15 of the FCC Rules and with Industry Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation. Changes or modifications made to this equipment not expressly approved by Withings may void the FCC authorization to operate this equipment.

Radiofrequency radiation exposure information:

The radiated output power of the device is far below the FCC radio frequency exposure limits. Nevertheless, the device shall be used in such a manner that the potential for human contact during normal operation is minimized.

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EN I This device complies with the following normative documents: COUNCIL DIRECTIVE 93/42/EEC of 14 June 1993 concerning medical devices as amended by Directive 2007/47/EC

EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes. Reference to standards contd. | EN ISO 14971: Medical devices — Application of risk management to medical devices. | IEC/EN 60601-1: Medical electrical equipment — Part 1: General requirements for basic safety and essential performance. | EN 1060-3: Non-invasive sphygmomanometers, Part 3: Supplementary requirements for electromechanical blood pressure measuring systems. | EN 1060-4: Non-invasive sphygmomanometers. Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers. | IEC/EN 60601-1-11: General requirements for basic safety and essential performance — Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. | IEC 80601-2-30: Medical electrical equipment — Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers. | EN 300 328: Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2.4 GHz ISM band and using wide band modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the Directive (2014/53/EU). | EN 301 489-1: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment

and services; Part 1: Common technical requirements. | EN 301 489-1: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements | IEC/EN60601-1-2: Medical electrical equipment: Part 1-2: General requirements for basic safety and essential performance-collateral standard electromagnetic compatibility | EN 55011: Industrial, scientific and medical equipment — Radio-frequency disturbance characteristics — Limits and methods of measurement. | FCC part B 15B: Electromagnetic Compatibility. | FCC Rule Part: 15.247 Cat: DSS (Bluetooth) FCC Rule Part: 15.247 Cat: DTS (BT4.0). | EN ISO 10993-1: Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process. | EN ISO 10993-5: Biological evaluation of medical devices — Part 5: Tests for in vitro cytotoxicity | EN ISO 10993-10: Biological evaluation of medical devices. Tests for irritation and skin sensitization | IEC/EN 60601-1-6: Medical electrical equipment — Part 1-6: General requirements for basic safety and essential performance — Collateral standard: Usability | ANS/AAMI/ISO 81060-2 Non-invasive sphygmomanometers Part 2: Clinical validation of automated measurement type | EN ISO 15223-1: Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements | IEC 60601-2-47:2012: Medical electrical equipment — Basic safety and essential performance of ambulatory electrocardiograph systems | IEC 57:2012: Testing and reporting performance results of cardiac rhythm and ST segment measurement algorithms

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EN ISO 13485: Medical devices – Quality management systems – Requirements for regulatory purposes. Reference to standards contd. | EN ISO 14971: Medical devices – Application of risk management to medical devices. | IEC/EN 60601-1: Medical electrical equipment - Part 1: General requirements for basic safety and essential performance. | EN 1060-3: Non-invasive sphygmomanometers, Part 3: Supplementary requirements for electromechanical blood pressure measuring systems. | EN 1060-4: Non-invasive sphygmomanometers. Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers. | IEC/EN 60601-1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment. | IEC 80601-2-30: Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers. | EN 300 328: Electromagnetic compatibility and Radio spectrum Matters (ERM); Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz ISM band and using wide band

modulation techniques; Harmonized EN covering the essential requirements of article 3.2 of the Directive (2014/53/EU). | EN 301 489-1: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements. | EN 301 489-17: Electromagnetic compatibility and Radio spectrum Matters (ERM); ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements | IEC/EN60601-1-2: Medical electrical equipment: Part 1-2: General requirements for basic safety and essential performance-collateral standard electromagnetic compatibility | EN 55011: Industrial, scientific and medical equipment – Radio-frequency disturbance characteristics – Limits and methods of measurement. | FCC part B 15B: Electromagnetic Compatibility. | FCC Rule Part: 15.247 Cat: DSS (Bluetooth) FCC Rule Part: 15.247 Cat: DTS (BT4.0). | EN ISO 10993-1: Biological evaluation of medical devices -- Part 1: Evaluation and testing within a risk management process. | EN ISO 10993-5: Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity | EN ISO 10993-10: Biological evaluation of medical devices. Tests for irritation and skin sensitization | IEC/EN 60601-1-6: Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability | ANSI/AAMI/ISO 81060-2 Non-invasive sphygmomanometers Part 2: Clinical validation of automated measurement type | EN ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements | IEC 60601-2-47:2012: Medical electrical equipment - Basic safety and essential performance of ambulatory electrocardiograph systems | IEC 57:2012: Testing and reporting performance results of cardia rhythm and STsegment measurement algorithms

Manufacturer/Fabricant/Hersteller

WITHINGS,
2 rue Maurice Hartmann,
92130 Issy-les-Moulineaux, FRANCE
+33141460460

1282

IP22

Ingress of water or particulate matter

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Complies with waste electrical and electronic equipment directive

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Frequency band :
2402 - 2480 MHz 12822460
Maximum output power : 4.45 dBm

EN Withings hereby declares that the device Withings BPM Core is in conformity with the essential requirements and other relevant requirements of Directive 2014/53/EU. The declaration of conformity can be found at: withings.com/compliance

UK Responsible Person:
Emergo Consulting (UK) Ltd c/o Cr360 - UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ United Kingdom

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AUS Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Australian sponsor:
Emergo Australia
Level 20, Tower II, Darling Park, 201 Sussex Street, Sydney, NSW 2000 Australia

Need help?

Besoin d'aide ? Brauchen Sie Hilfe?

➔ withings.com/support

AW_WPM04_QIG_INTER_A

WITHINGS

Withings BPM Core | Smart blood pressure monitor with ECG & digital stethoscope | V1.2 — May. 2022

© 2022 Withings. All rights reserved.