## **Scan Monitor Instructions for Use**

### **Withings**

2 rue Maurice Hartmann  
92130 Issy-les-Moulineaux, FRANCE  
[Withings.com](http://Withings.com)

#### **1. INDICATIONS FOR USE**

Scan Monitor is a medical device composed of software and the dedicated hardware of a reusable wrist applied device, ScanWatch, which incorporates a lead I electrocardiograph (ECG).

Scan Monitor is intended to measure heart rate and to detect atrial fibrillation with a single channel electrocardiogram. It is intended for intermittent measurements.

Scan Monitor measures, transfers, records, and displays a single channel electrocardiogram similar to a lead I ECG. Scan Monitor's software detects the presence of atrial fibrillation (AFib) or sinus rhythm on a classifiable ECG waveform, calculates the heart rate, the P wave duration, the PR interval, the QT interval, the QT corrected interval, and the QRS duration.

Scan Monitor is indicated for use in adults (18 years of age or older).

#### **2. USING THE SCAN MONITOR**

##### **a. Set-Up**


- Withings ScanWatch is a connected watch, distributed without ECG function
- Scan Monitor is activated during the Association of ScanWatch and the Withings app.
- Scan Monitor is only available for Withings ScanWatch, paired with a smartphone with iOS 15 or later, or Android 9 or later.

##### **b. On-boarding**

- Open the Withings app.
- In the timeline tab, tap “+”, then select “Electrocardiogram (ECG)”
- Follow the onscreen instructions.

## Scan Monitor

- You may exit on-boarding at any time by pressing “Cancel”.
- c. How to take an ECG recording
- Make sure your ScanWatch is snug on the wrist you selected. Tap Devices > More settings> Device Orientation.
- Rest your arms on a table and hold the top electrode with your second hand (with your thumb and index for example, as shown on fig. 1). You do not need to press the bezel during the session.
- Make sure to follow the gesture indicated in “the best practice” article within the mobile application in order to obtain an ECG with a good signal quality

*Fig. 1 - How to place second hand to get an ECG measurement*

- The recording starts after the first vibration.
- The recording lasts for 30 seconds.
- End of measurement is confirmed by a double vibration.

## 3. ECG Analysis

- After a successful reading, you will receive one of the following classifications on your Scan Monitor:
  - Sinus Rhythm: a sinus rhythm result means the heart is beating in a uniform pattern between 50-100 BPM.
  - Atrial Fibrillation: an AFib result means the heart is beating in an irregular pattern above 50 BPM.
  - Inconclusive: An inconclusive result means the recording can't be classified. This can happen for many reasons such as not resting your arms on a table during a recording, or your ScanWatch is too loose. Certain physiological conditions may prevent a small percentage of users from creating enough signal to produce a good recording.

## Scan Monitor

- Low or High Heart Rate: A heart rate under 50 BPM or over 100 BPM affects the ECG app's ability to check for AFib and the recording is considered inconclusive.
- After an ECG recording is complete, the ECG data is analyzed to determine if it is 30-second long, and, if so, if either Sinus Rhythm or AFib is present, or if an Inconclusive result is warranted.
- The ECG recording results on a detailed display of the result on the ScanWatch and on the Withings app.
- The classification of the ECG recording is for information use only. It is meant to supplement, but not replace traditional methods of diagnosis. If you are experiencing any symptoms or have concerns, contact your physician. If you believe you are experiencing a medical emergency, you should contact emergency services.
- A result of Sinus Rhythm means your heart rate is between 50 and 100 beats per minute (bpm) and is beating in a regular pattern.
- Inconclusive ECG results may be caused by too many artefacts, by movements (signal of poor quality), by too much noise parasitizing due to the proximity with an electrical device generating strong electromagnetic fields or you may have an arrhythmias other than AFib or your heart is below 50bpm or above 100bpm. A small percentage of people may have certain physiological conditions preventing the user from creating enough signal to produce a good recording. You can learn more about Inconclusive ECG results during on-boarding, by accessing educational information in the ECG area of the Withings app on your smartphone.
  - A heart rate can be low because of certain medications or if electrical signals are not properly conducted through the heart. Training to be an elite athlete can also lead to a low heart rate.
  - A heart rate can be high because of exercise, stress, nervousness, alcohol dehydration, infection, AFib, or another arrhythmia.
  - If you receive an Inconclusive result due to a poor recording, you might try to re-record your ECG. You can review how to take an ECG during on-boarding or by tapping on "Take a Recording" in the ECG area of the Withings app on your smartphone.
- All ECGs are synced to the Withings app on your Android or iOS smartphone. You may use the Withings app to share your ECG with a clinician.

## 4. TECHNICAL SPECIFICATIONS

Product Name	Scan Monitor
ECG Sensor	1 stainless steel electrode (back case electrode) & 1 Titanium Electrode (bezel electrode) or 2 stainless steel electrodes coated with TiCN/Au
Operating conditions	+5°C to 35°C 15 to 90% Relative humidity 700-1060hPa
Storage and transport conditions	-25°C to 70°C

## Scan Monitor

	20 to 90%
Power Source	3.6 Vdc Lithium ion battery
Wireless transmission	Wi-Fi and BLE
Measurement range (Heart Rate)	30bpm to 230bpm
Display	LED on the watch and in-app

## 5. SAFETY AND PERFORMANCE

The Scan Monitor software's ability to accurately classify an ECG recording into Afib and sinus rhythm was tested against the IEC 60601-2-47\* standard and clinical data of 168 subjects.

The Scan Monitor demonstrated 98.16 % sensitivity in classifying AFib and 97.20 % specificity in classifying sinus rhythm in classifiable recordings of 5 public databases.

Rhythm classification from Scan Monitor was compared to ECG labeled by cardiologist on the clinical data. During this clinical study, 20.2% of recordings were inconclusive and not classifiable as either sinus rhythm or AFib. The Scan Monitor reached 100.0 % sensitivity in classifying AFib and 100.0 % specificity in classifying sinus rhythm among the remaining classifiable recordings. These results reflect use in a controlled environment. Real world use of the Scan Monitor may result in a greater number of strips being deemed inconclusive.

The waveform output from the software was tested in accordance with IEC 60601-2-47\* standard by a beat-to-beat QRS detection. This detection reaches a F1-score of at least 99,19% on all the datasets with the exception of NSTDB (90,65% because of the digitally added noise).

*\*\*IEC 60601-2-47:2012 : Requirements for the Basic Safety and Essential Performance of Ambulatory Electrocardiographic Systems.*

### 5.1. RF Statement

Guidance and manufacturer's declaration-electromagnetic emissions		
The Scan Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Scan Monitor should assure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment-guidance
CE emissions CISPR11	Group 1	The Scan Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RE emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	The Scan Monitor is suitable for use in all establishments, including domestic
Voltage fluctuations/ Flicker emissions IEC 61000-3-3		establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

### 5.2. Declaration- electromagnetic emissions and immunity for equipment and systems that are not life-supporting and are specified for use only in a shielded location

## Scan Monitor

Declaration- Electromagnetic immunity			
The Scan Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Scan Monitor should assure that it is used in such an environment			
Immunity test	IEC 60601 Test Level	Compliance level	Electromagnetic environment-guidance
Conducted	3 Vrms	N/A	N/A
RF IEC 61000-4-6	150kHz to 80MHz	N/A	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the equipment or system including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol.
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8kV Air: ±2kV, ±4kV, ± 8kV, ±15kV	Contact: ± 8kV Air: ±2kV, ±4kV, ± 8kV, ±15kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/ burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	N/A	A mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	-5% UT (95% dip in UT) for 0.5 cycles, -40% UT (60% dip in UT) for 5 cycles, -70% UT (30% dip in UT) for 25 cycles, -5% UT (95% dip in UT) for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power main interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/ 60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment

## 6. TROUBLESHOOTING

If you experience difficulties in operating your Scan Monitor, refer to the troubleshooting guide below.

### 6.1 Troubleshooting for ECG measurements

Problem 6.1.1: I cannot get the Scan Monitor to take an ECG reading.

Solutions:

- Make sure your wrist and your ScanWatch are clean and dry. Water and sweat can cause a poor recording.
- Ensure that your ScanWatch, arms, and hands remain still during recordings.
- Ensure that you have completed all of the on-boarding steps in the Withings app on your smartphone.

Problem 6.1.2: I have an inconclusive measure. It looks like the ECG measure presents a lot of artefacts, noise, or interference in my recording.

Solution:

- Rest your arms on a table while you take a recording. Try to relax and don't move too much.
- Tighten the band so that the back of the watch is in contact with the skin of the wrist. When moving the watch slightly, the skin in contact with it should move with it.
- Move away from any electronics that are plugged into an outlet to avoid electrical interference.

Problem 6.1.3: The ECG waveforms appear upside down.

Solution:

- The device orientation may be set to the wrong wrist. On your smartphone, go to the Withings app. Tap Devices > More Settings > Device Orientation.

All data recorded during an ECG measurement is saved to the Withings app on your smartphone. If you choose to, you can share that information by creating a PDF.

## 7. Cleaning and disposal

### 7.1 Cleaning ScanWatch

- Use a lint-free cloth moistened with warm water to clean the top housing and casing of your ScanWatch
- Run the wrist band under warm water and rub it with hypoallergenic soap to clean it
- Dry the wrist band with a soft cloth

## 7.2 Disposal

The ScanWatch is classified as an electrical and electronic equipment. Such items should not be mixed with general household waste and must be taken to dedicated collection points at the end of their working life for proper treatment, recovery and recycling.

## 8. Contraindications

Scan Monitor is not intended to continuously monitor vital signs in critical conditions or where the nature of the variations is such that it could result in immediate danger to the patient. Scan Monitor does not provide alarms.

Scan Monitor is not intended to provide or to replace a diagnosis by a physician or qualified health care professionals. Vital signs measurements, such as those taken with this device, cannot identify all diseases.

Scan Monitor is not indicated for use:

- On patients with a pacemaker or other implanted electronic device
- On patients diagnosed with an arrhythmia other than atrial fibrillation
- For out-of-hospital transport by emergency medical services
- During medical procedure

## 9. Warnings

Regardless of the measurement taken using this device, you should consult your doctor immediately if you experience symptoms that could indicate acute disease.

Scan Monitor cannot detect heart attacks or ischemic heart conditions. If you ever experience chest pain, pressure, tightness or what you think is a heart attack, call emergency services.

DO NOT self-diagnose or self-medicate on the basis of this device without consulting your doctor. In particular, do not start taking any new medication or change the type and/or dosage of any existing medication without prior approval.

Misapplication of the ScanWatch with excessive pressure and for a prolonged period can cause pressure injury.

Avoid application of the ScanWatch to a wrist with poor skin integrity.

DO NOT sterilize by irradiation, steam or ethylene oxide- refer to cleaning and disinfection instructions. Use of cleaning agents other than specified may damage the device.

DO NOT use the ScanWatch if it is damaged- use of a damaged device could cause patient injury or equipment failure.

## Scan Monitor

DO NOT take recordings when ScanWatch is in close vicinity to strong electromagnetic fields (e.g. electromagnetic anti-theft systems, metal detectors).

DO NOT take recordings during a medical procedure (e.g., magnetic resonance imaging, diathermy, lithotripsy, cautery and external defibrillation procedures).

DO NOT take recordings when ScanWatch is outside of the operational temperature range (5°C - 35°C) indicated in the ScanWatch user manual and humidity range of 15% to 90% relative humidity.

DO NOT use to diagnose heart-related conditions.

DO NOT use with a cardiac pacemaker, ICDs, or other implanted electronic devices.

DO NOT take a recording during physical activity.

### NOTE:

The device is conforming to IP22 requirements.

**CAUTION:** Withings does not guarantee that you are not experiencing an arrhythmia or other health conditions when the Scan Monitor labels an ECG as Sinus Rhythm. You should notify your physician if you detect possible changes in your health.

## 10. EQUIPMENT SYMBOLS

Symbol	Description
	Manufacturer
	Serial Number
	Temperature limit
	Humidity limitation
	Atmospheric Pressure limitation

Scan Monitor

<b>UDI</b>	Unique Device Identifier
	Consult instructions for use
	Follow instructions for use
	Importer
<b>CH REP</b>	Swiss Authorised Representative

**Australian Sponsor:**

Emergo Australia, Level 20, Tower II  
Darling Park, 201 Sussex Street- Sydney, NSW 2000  
Australia

**UK Responsible Person:**

Emergo Consulting (UK) Limited  
c/o Cr360 – UL International, Compass House,  
Vision Park Histon, Cambridge CB24 9BZ, Royaume-Uni

MedEnvoy Switzerland

Gotthardstrasse 28

6302 Zug

Switzerland

## **Scan Monitor Mode d'emploi**

Withings  
2 rue Maurice Hartmann  
92130 Issy-les-Moulineaux, FRANCE  
[Withings.com](http://Withings.com)

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