EN
FR
DE

WITHINGS SCANWATCH NOVA

PRODUCT GUIDE

USER GUIDE

EN - NON-US REGIONS

This leaflet explains how to use the Withings ScanWatch NOVA.

WITHINGS, 2 rue Maurice Hartmann,
92130 Issy-les-Moulineaux, FRANCE +33141460460

Important notice

Before using Withings ScanWatch NOVA, review the information in this guide. You may also find this guide online at: <https://www.withings.com/guides>.

Retain this documentation for future reference. Installation instructions are available with the Quick Start Guide provided with this Product Guide.

Please contact Withings if you need assistance, setting up, using or maintaining the device or need to report unexpected operations or events. Any serious incident that has occurred when using Withings ScanWatch NOVA should be reported to Withings and appropriate authorities in your country of residence.

Disclaimer

Information in this guide may change without notice.

Intended Use (Non-US Regions)

Withings ScanWatch NOVA is a medical device composed of software and the dedicated hardware of a reusable wrist applied device, which incorporates a single-lead electrocardiograph (ECG) and a reflectance photoplethysmograph (PPG). It is intended for intermittent measurements.

Withings ScanWatch NOVA measures, transfers, records, and displays lead I of an ECG. It calculates the heart rate and detects the presence of AF or sinus rhythm on a classifiable ECG

waveform.

Withings ScanWatch NOVA measures, transfers, records and displays pulse rate to identify episodes of irregular heart rhythm suggestive of AF and provides a notification to the user. It can be used to supplement a clinician's decision to screen for possible AF.

Withings ScanWatch NOVA is indicated for over-the-counter (OTC) use in adults.

Withings ScanWatch NOVA is not recommended for users with other known arrhythmias.

Contraindications

DO NOT use the ECG and irregular heart rhythm notification features if you have a cardiac pacemaker, implantable cardiac-defibrillator, or other implanted electronic device.

This irregular heart rhythm notification feature cannot detect heart attacks. If you ever experience chest pain, pressure, tightness, or what you think is a heart attack, call emergency services.

The Irregular Rhythm Notification Feature is not constantly looking for AF and should not be relied on as a continuous monitor. This means the feature cannot detect all instances of AF, and people with AF may not get a notification.

This assessment carried out by this product is not a diagnosis.

The data displayed by the Withings ScanWatch NOVA is intended for informational use only, and is not intended to replace traditional methods of diagnosis or treatment.

Regardless of the measurements taken using this device, you should consult your doctor immediately if you experience symptoms that could indicate you are experiencing a sudden and/or severe change in health.

DO NOT self-diagnose or self-medicate on the basis of this device without consulting your doctor. In particular, **DO NOT** start taking any new medication or change the type and/or dosage of any existing medication without prior approval from a healthcare professional.

The device is not intended to continuously monitor vital physiological parameters in critical conditions or where the nature of the variations is such that it could result in immediate danger to the user.

The device does not provide alarms.

Warnings

Before using this equipment, be sure to read this manual carefully.

Your device uses Bluetooth technology. Keep your Withings ScanWatch NOVA device more than 15 cm away from any implanted medical devices, and other devices that may be affected. Stop using your Withings ScanWatch NOVA if you suspect any interference with your medical device, and consult your physician or your medical device manufacturer as needed.

A number of factors can affect the accuracy of the pulse rate data used to identify episodes of irregular heart rhythm suggestive of AF. These include factors like body motion, hand and finger movements, wrist tattoos and decreased blood flow to your wrist (which can be caused by cold temperatures). Exposing the device to prolonged lint or dust may reduce its lifetime or damage it.

Inspect electrodes for warping, surface damage or corrosion and check for any other form of damage. Damaged electrodes or sensors may lead to incorrect measurements. **DO NOT** use Withings ScanWatch NOVA if it is damaged. Use of a damaged device could cause injury or equipment failure. **DO NOT** try to repair or modify this equipment yourself. **DO NOT** open or disassemble the device for battery replacement.

Use of cleaning agents other than specified in this manual may damage the device. **DO NOT** use accessories, detachable parts, or materials not described in this manual. Use of parts and components other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. **DO NOT** interconnect this equipment with other equipment not described in this manual. **DO NOT** take recordings when Withings ScanWatch NOVA is in close vicinity to strong electromagnetic fields (e.g., electromagnetic anti-theft systems, metal detectors). **DO NOT** take recordings during a medical procedure (e.g., magnetic resonance imaging, diathermy, lithotripsy, cautery and external defibrillation procedures). **DO NOT** take recordings when Withings ScanWatch NOVA is outside of the operational temperature and relative humidity ranges indicated in this manual. Otherwise, measurement results may be impacted.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise this could result in degradation of performance of the Withings ScanWatch NOVA.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

The power cord of the charger may cause strangulation due to excessive length. Keep it away from children and pets.

The battery inside the watch will stop charging when the temperature is less than -10°C (+/-5°C or) or over 50°C (+/-5°C).

If storage conditions are different from the operating conditions indicated in the technical specifications, wait 30 minutes before using the device.

Accessory

The Withings ScanWatch NOVA is to be used with a software accessory, the Withings ScanWatch 2 companion app. It allows the installation of the watch, activates its medical features, and displays the results of performed measurements. It is part of the Withings App.

Getting started


Box Content

ScanWatch NOVA

Quick Start Guide
and Product Guide

Charger and cable

Fluoroelastomer
wristband


Hammer

Pin punch


Holder


Note: To recharge the battery, use the provided power cord.

Band toolkit

Note: This toolkit is only suitable for a metal wristband with a folding clasp.


Adjust the wristband on your wrist to estimate the number of mesh links to be removed. Note that one of the links is smaller in order to adjust the bracelet more finely.


Position the wristband in the HOLDER. the arrows inside the wristband should point downward.


Band toolkit

Use the pin punch and the hammer to push out the pin that supports the mesh you want to remove.

Remove the second pin. Remove the mesh. Repeat if necessary.

Reassemble the wristband by replacing one of the pins, then the other (in the opposite direction of the arrows) with the hammer.

How to use Withings ScanWatch NOVA

Set-Up

- Follow the information in the Quick Start Guide and in the Withings App to set-up your device.
- To wake up your device for the first time, please press the crown. If your product is not starting, consider charging it for a minimum of two hours and try again.
- Once your product is installed and configured, it will appear in your devices list in settings.
- You may have to follow some tutorials in the Withings App in order to activate some measurements. The ECG and Irregular Heart Rhythm Notification features require you to follow onboarding tutorials for their activation.

What is the Irregular Heart Rhythm Notification feature?

To activate this feature, you must follow the onboarding tutorial in the Withings App. Withings ScanWatch NOVA can look for signs of atrial fibrillation. Every 10 minutes, your watch looks for irregular patterns in your heart rhythm in the background. Once the feature is activated, an irregular heartbeat scan is launched every 10 minutes without any input from your side. As soon as an irregular heart rhythm suggestive of atrial fibrillation is detected, you'll be notified on your device and on your phone. In that case, your watch will suggest taking an ECG recording.


Note: this text ("Signs of Afib") will scroll across the watch screen.

♥ Signs of AFib
Detected
 ⚠ 3 days that month

How to take an ECG recording

- Make sure your Withings ScanWatch NOVA is snug on the wrist you selected during the tutorial. To see which wrist you have chosen, on the Withings App tap Devices > All settings > ECG.
- Rest your arms on a table and hold the top electrode (bezel) with your second hand (with your thumb and index for example), as shown on the figure thereafter. You do not need to press the bezel during the session.
- Make sure to follow the gesture indicated in "the best practice" article within the Withings App in order to obtain an ECG with a good signal quality.
- The recording starts after the first vibration.
- The recording lasts for 30 seconds.
- End of measurement is confirmed by the second vibration.


ECG Outputs

ECG classification:

After an ECG recording, you will see one of the following classifications for the recording on the watch screen (cf. left picture of the figure below) and in the Withings App (cf. right picture of the figure below):


- Low Heart Rate (heart rate < 50 bpm)
- High Heart Rate (heart rate > 150 bpm)
- Sinus Rhythm (heart rate between 50-99 bpm)
- High Heart Rate (No signs of AFib) (heart rate between 100-150 bpm)
- Atrial Fibrillation (heart rate between 50-99 bpm)
- Atrial Fibrillation — High Heart Rate (heart rate between 100-150 bpm)
- Inconclusive
- Poor Recording

After the ECG recording you will also see your average heart rate and your ECG tracing.


Low Heart Rate

A Low Heart Rate result means your heart is beating less than 50 beats per minute (bpm). This recording cannot be classified by the device. A low heart rate can happen if electrical signals are not properly conducted through the heart. Some medicines can also cause a low heart rate. Talk to your doctor if you have questions about your ECG recording.

High Heart Rate

A High Heart Rate result means your heart is beating above 150 beats per minute (bpm). This recording cannot be classified by the device. Many different things can cause a high heart rate. A heart rate may be high because of exercise, stress, dehydration, infection, AFib, another arrhythmia or another cause. Talk to your doctor if you have questions about your ECG recording.

Sinus Rhythm

A Sinus Rhythm result means your heart rate is between 50 and 99 beats per minute (bpm) and is beating regularly.

High HR


High Heart Rate (No signs of AFib)

A High Heart Rate (No signs of AFib) result means the heart rate is beating between 100 and 150 beats per minute (bpm) and does not show any signs of Atrial Fibrillation.

Many different things can cause a high heart rate. A heart rate may be high because of exercise, stress, dehydration, infection, an arrhythmia, or another cause. Talk to your doctor if you have questions about your ECG recording.

Atrial fibrilla

Note: this text (“Atrial Fibrillation”) will scroll across the watch screen.

Atrial Fibrillation

An Atrial Fibrillation result means the heart rate is between 50 and 99 beats per minute (bpm) and is beating irregularly. If you have not been diagnosed with AFib before, you should talk to your doctor.

Atrial fibrilla

Note: this text (“Atrial Fibrillation”) will scroll across the watch screen.

Atrial Fibrillation - High Heart Rate

An Atrial Fibrillation — High HR result means your heart rate is beating between 100 and 150 beats per minute and is beating with an irregular pattern. If you have not been diagnosed with AFib before, you should talk to your doctor.

Inconclusive

Inconclusive

An Inconclusive result means that the signal cannot be classified as Sinus Rhythm or Atrial Fibrillation, although the quality of the recording is good. This may be due to various conditions, including but not limited to other arrhythmias or other heart conditions. Talk to your doctor if you have questions about your ECG recording.

Poor Recording

A Poor Recording result means the recording quality is low and the ECG cannot be classified. Some things that can cause this type of result are: excessive movement that causes a signal of poor quality; OR proximity to an electrical device that generates strong electromagnetic fields; OR not following best practices for gestures to be adopted ; A small percentage of people may have certain physiological conditions that prevent them from creating enough signal to produce a quality recording. You may try to re-record your ECG. You can review how to take an ECG during setup or by tapping Take a Recording in the ECG section of the Withings app on your smartphone. If you think you may be having a heart attack (myocardial infarction) or are facing a medical emergency, call the emergency services.

- The classification of the ECG recording is for information use only. It is meant to supplement, but not replace traditional methods of diagnosis. If you are experiencing any symptoms or have concerns, contact your physician. If you believe you are experiencing a medical emergency, you should contact emergency services.
- All ECGs are synced to the Withings app on your Android or iOS smartphone. You may use the Withings app to share your ECG with a clinician.

Sharing ECG Results

How to share your ECG with your doctor:

Once the ECG measurement feature is activated (see the section “How to use Withings ScanWatch NOVA” of this product guide), the results of the ECG can be shared with a physician as a PDF, generated by the Withings app.

The pdf includes the following information:

- The ECG strip for the lead I.
- The average heart rate
- The classification of the ECG according to the detection algorithm.

Troubleshooting

If you cannot fix the problem using the following troubleshooting instructions, please contact Withings or go to: withings.com/support.

Troubleshooting for ECG measurements

Problem	Solution
I cannot get the Withings ScanWatch NOVA to take an ECG reading.	- Make sure your wrist and your Withings ScanWatch NOVA are clean and dry. Water and sweat can cause a poor recording. - Ensure that your Withings ScanWatch NOVA, arms, and hands remain still during recordings. - Ensure that you have completed all of the on-boarding steps in the Withings app on your smartphone.
I have an inconclusive measure. It looks like the ECG measure presents a lot of artifacts, noise, or interference in my recording.	Rest your arms on a table while you take a recording. Try to relax and don't move too much. - Tighten the band so that the back of the watch is in contact with the skin of the wrist. When moving the watch slightly, the skin in contact with it should move with it. - Move away from any electronics that are plugged into an outlet to avoid electrical interference.
The ECG waveforms appear upside down.	The device orientation may be set to the wrong wrist. On your smartphone, go to the Withings app. Tap Devices > More Settings > Device Orientation. All data recorded during an ECG measurement is saved to the Withings app on your smartphone. If you choose to, you can share that Information by creating a PDF.

Technical specifications

Product Name: Withings ScanWatch NOVA
Model: HWA10
ECG Sensor: For hwa10 38 mm and hwa10 42 mm: 1 stainless steel electrode (back case electrode) & 1 Titanium Electrode (bezel electrode)
For hwa10 Rosegold: 2 stainless steel electrodes coated with TiCN/Au or 2 stainless steel electrodes for hwa10 Nova 42 mm.
PPG Sensor: Red light (wavelength: 660nm) Infrared light (wavelengths: 850nm and 940nm) Green (wavelength: 530nm)
Accelerometer sensitivity: +/-4g
Battery Operated: Up to 30 days on a single charge.
Power Source: Powered internally by a 3.87 Vdc Lithium ion battery (use the charging dock and USB-C charging cable included) and a DC 5V power adapter with a minimum power of 3W and a maximum power of 5W (not included).
Heart Rate Measurement Range (ECG feature): 30 bpm to 220 bpm
Heart Rate (Accuracy): Bias: -0.97 bpm, standard deviation: 3.3 bpm
Display: OLED screen on the watch
Weight Withings ScanWatch NOVA: 63 g (42 mm)
Operation time: 30 days typical use on a single charge
Parts in contact with the skin: Fluoroelastomer (wristband), mineral glass (lens), stainless steel (case), stainless steel (coated with TiCN/Au (case), titanium (case), sap-phire glass (front glass), thermoplastic polyester elastomer (compression rings on back case).
Transport and storage conditions: Temperature range: -25°C to 70°C, Humidity range: 20 to 90%, Pressure range: 700-1060hPa
Operation conditions: Temperature range: 5°C to 35°C, Humidity range: 15% - 90%, Pressure range: 700-1060hPa
If storage conditions are different from these operating conditions, wait 30 minutes before using the device.
Wireless Transmission: BLE
IP Protection level: IP22
Mode of operation: Continuous Operation
Expected product life: 3 years
Type BF-Applied Part: back case

Wireless information

Wireless Specifications

Wireless Technology	Bluetooth BLE
Version	Supported BT5.1
Operation Frequency	2402MHz - 2480MHz
Transmission Power	+8dBm (max)
Modulation	GFSK
Receiver sensitivity	-96dBm

The wireless communication of the Withings ScanWatch NOVA is supported by a BLE communication. This communication is established between the Withings Scan-Watch NOVA and the Withings App. The communication between the Withings ScanWatch NOVA and the Withings App is encrypted through an exchange of a paired key.

The communication latency between the Withings ScanWatch NOVA and the Withings App is inferior to 10 seconds. A maximum operating distance of 5 meters allows a latency inferior to 10 seconds. The communication security is implemented by default (encrypted communication). The Withings App shall be downloaded from official stores (App store and Google Play store) and smartphones shall be up to date. Versions supported by the Withings App are iOS 15 or later or Android 9 or later. In case of communication failure, you should follow the related troubleshooting. Measurements are stored within the Withings ScanWatch NOVA. The communication between the Withings ScanWatch NOVA and the Withings App is not modified with sources of interference signals located within 5 meters.

This wireless coexistence has been tested in accordance with the following standards:
ANSI C63.27:2021 and AAMI TIR69:2017 (R2020)

Electromagnetic disturbances have been tested in accordance with the standard IEC 60601-1-2:2014/A1:2020

Cleaning & disposal

Cleaning of Withings ScanWatch NOVA

- Use a lint-free cloth moistened with warm water to clean the glass and casing of your Withings ScanWatch NOVA.
- Run the wristband under warm water and rub it with hypoallergenic soap to clean it
- Dry the wrist band with a soft cloth

Notes:

- Avoid exposure to any chemical products such as solvents, detergents, perfumes, or cosmetics. They may damage the wristband or the casing of the watch.
- If the side button of your Withings ScanWatch NOVA appears to be stuck, we advise you to rinse the casing of your watch with warm water and dry it using a lint-free cloth.
- If the altimeter opening of Withings ScanWatch NOVA becomes dirty or clogged with debris or sweat, it can prevent the watch from accurately measuring elevation. If this occurs, please clean the case of ScanWatch NOVA with warm, soapy water and then allow it to dry. The watch should also be cleaned anytime that you work out with it on
- Withings does not recommend using an ultrasonic cleaner for Withings ScanWatch NOVA.

Avoid application of the Withings ScanWatch NOVA on a wrist with poor skin integrity.

In case of skin irritation, we advise you to contact a dermatologist before resuming use of your device.

Skin reactions may be caused by:

- An allergic reaction to the materials used in the back casing of the watch,
- A photosensitizing medication.

Note that this device is compliant to IEC 62471 related to the photobiological safety of lamps and lamp systems.

European Union - Declaration of conformity

Withings hereby declares that the Withings ScanWatch Nova and Withings ScanWatch 2 companion app are in conformity with the general safety and performance requirements of the Medical Devices Regulation (EU) 2017/745. The full text of the EU declaration of conformity can be found at: withings.com/compliance

Safety & performance

Performance and intended clinical benefit

The Withings ScanWatch NOVA contributes to an early detection of atrial fibrillation, facilitates patient screening and helps physicians treat the disease earlier and thus attenuates the impact of the disease on the population. It also contributes to an increased awareness of physiological conditions.

Withings ScanWatch NOVA's ability to accurately classify an ECG recording into Atrial Fibrillation (AF) and Sinus Rhythm (SR) categories was tested in a clinical study including 626 patients. Rhythm classification by Withings ScanWatch NOVA was compared to the review by cardiologists of simultaneously recorded 12-lead ECG. 89.1% of recordings were conclusive (i.e. neither "inconclusive" nor "poor recording", and with a heart rate between 50 and 150 bpm). On conclusive recordings, the sensitivity in classifying Atrial Fibrillation was 99.7% (lower bound of the 95% confidence interval: 98.1%) and the specificity in detecting normal sinus rhythm was 99.8% (lower bound of the 95% confidence interval: 99.1%). In addition, ScanWatch NOVA accurately classifies AF and SR into subgroups by heart rate < or ≥ 100 bpm.

The mean difference between the algorithm calculated heart rate from Withings ScanWatch NOVA's strips (median over 30 seconds) and the heart rate measured by cardiologists on the 12-lead ECGs was -0.97 bpm (standard deviation: 3.3 bpm).

The quality of the ECG waveforms was assessed by the visibility and polarity of the P, QRS and T waves. The accuracy of the 6 endpoints ranged between 97.8% and 100%. In addition, the mean difference (standard deviation) for PR, QRS, and QT intervals was 3.4 (17.2) ms, -10.5 (15.9) ms, and 5.3 (23.8) ms respectively. Cardiologists also reviewed all the ECG strips generated via the Withings ScanWatch 2 companion app. They found that 92.6% of the recordings were

conclusive, and reached a sensitivity of 100% and a specificity of 97.9%.

The performance of the ECG function passed the applicable requirements on ECG waveform and heart rate measurements specified in ANSI/AAMI/IEC 60601-2-47:2012 Requirements for the Basic Safety and Essential Performance of Ambulatory Electrocardiographic Systems and ANSI/AAMI EC57:2012 (R2020) Testing And Reporting Performance Results Of Cardiac Rhythm And ST Segment Measurement Algorithms.

The ability of Withings ScanWatch NOVA to accurately classify a PPG recording into Atrial Fibrillation and Sinus Rhythm categories was assessed through a clinical study including 166 patients. Rhythm classification was compared to the review of simultaneously recorded 12-lead ECGs by trained cardiologists. On classifiable signals (i.e. neither "noisy" nor "other arrhythmia"), the positive predictive value for AF detection was 0.864.

RF statement

Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in the accompanying documents. Portable and mobile RF communications equipment can affect medical electrical equipment.

Safety and performance

Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment — guidance
CE emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RE emissions CISPR11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	

Declaration – electromagnetic immunity

Declaration- electromagnetic emissions and immunity for equipment and systems that are not life-supporting and are specified for use only in a shielded location

Withings ScanWatch NOVA is intended for use in the electromagnetic environment specified below. The customer or the user of the Withings ScanWatch NOVA should ensure that it is used in such an environment.

Immunity test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz	3 Vrms 150kHz to 80MHz	The EUT shall be placed on an insulating support 0.1 m above the ground reference plane (GRP). On all cables to be tested, coupling and decoupling devices shall be inserted. The coupling and decoupling devices shall be placed on the GRP, making direct contact with it at a distance of 0.1 m to 0.3 m from the EUT.
Radiated RF IEC 61000 4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the equipment or system including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Interference may occur in the vicinity of equipment marked with the following symbol.
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ±8 kV Air: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Contact: ±8 kV Air: ±2 kV, ±4 kV, ±8 kV, ±15 kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material the relative humidity should be at least 30%.

Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	N/A	The main power quality should be of the kind used in a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	line(s) to line(s): +/- 1kV	The main power quality should be of the kind used in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips: 0% residual for 0.5 cycle Voltage Dips: 0% residual for 1 cycle Voltage Dips: 70% residual for 25 cycles Voltage Interruptions: 0% residual for 250 cycles	N/A	The main power quality should be of the kind used in a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power main interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	The power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Federal Communication Commission interference statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

IMPORTANT NOTE:

Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 cm between the radiator & your body.

This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. Country Code selection feature to be disabled for products marketed to the US/CANADA

Disposal

Actuation of European directive 2012/19/EU, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal. The symbol applied on the device or its packaging means that at the end of its useful life, the product must not be disposed of with domestic waste. At the end of the device's useful life, the user must deliver it to a collection center for electric and electronic garbage, or return it to the retailer when purchasing a new device. Disposing of the product separately prevents possible negative environmental and health consequences deriving from inadequate disposal. It also allows the recovery of materials it is composed of to save energy and resources and avoid negative effects to the environment and health. In case of abusive disposal of devices by the user, will be applied administrative endorsements in compliance with current standards. The device and its parts must be disposed of as appropriate, in accordance with national or regional regulations.

Warranty

Withings™ Two (2) Years Limited Warranty - Withings ScanWatch NOVA

Withings, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, FRANCE (“Withings”) warrants the Withings-branded hardware product (“Withings Product”) against defects in materials and workmanship when used normally in accordance with Withings’ published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser (“Warranty Period”). Withings’ published guidelines include but are not limited to information contained in technical specifications, safety instructions or the quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error- free. Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product’s use.

AU: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

 Need help ?

➔ <https://support.withings.com/>

Security

Withings recommends that you add a passcode (personal identification number [PIN]), Face ID or Touch ID (fingerprint) to your phone to add a layer of security. It is important to secure your phone since you will be storing personal health information. Users should follow the password and two-factor authentication guidelines when logging into Withings App. For more information on password requirements refer to <https://support.withings.com/hc/articles/360012727737>. For more information on two-factor authentication refer to <https://support.withings.com/hc/en-us/articles/4409345943697> for iOS users and refer to <https://support.withings.com/hc/en-us/articles/21600361983633> for Android users. Users will receive email alerts in case of changes related to password, two-factor authentication and recovery code. Users will also receive additional software update notifications on the device via the Withings App, and updates are delivered wirelessly, encouraging rapid adoption of the latest security fixes. Users are able to see the currently installed firmware in the Withings App, under Devices > ScanWatch NOVA. This tab also indicates if an update is available.

Do not install the device on a smartphone that you do not own. Do not use a public Wi-Fi network you don't know. Use a trusted Wi-Fi network with your device. Use a secure channel when you are sharing ECG information with your doctor. Withings also recommends upgrading your Withings App when an upgrade is available. The Withings App is not intended to be used on a computer. No anti-virus software is needed. Only use official app stores to download the Withings application. In case of doubt, use the link [go.withings.com](https://support.withings.com/hc/en-us/articles/360010218558).

If needed, users can restore device configurations by following the factory reset procedure. For more information refer to <https://support.withings.com/hc/en-us/articles/360010218558>

Equipment Symbol Description

 Direct current

 Model Number

 Catalogue number

 Do not use if package is damaged and consult instructions for use

 Caution

 Type BF applied part (back case)

 Regulatory Compliance Mark (RCM) with ACMA's labelling requirements

 Follow instructions for use

 Medical device (Non-US users)

 Consult instructions for use or consult electronic instructions for use

 Manufacturer and date of manufacture

 Complies with waste electrical and electronic equipment directive 2012/19/EU

 Complies with the general safety and performance requirements of the Medical Device Regulation (EU) 2017/745

 Cardboard

 Indicates the authorised representative in Switzerland

 Indicates the importer in Switzerland

 Serial number

 Unique device identifier

 Temperature limit

 Atmospheric pressure limitation

 Humidity limitation

 Protected against solid foreign objects of 12,5 mm Ø and greater
Protection against vertically falling water drops when ENCLOSURE tilted up to 15°

 No charging device included with the radio equipment (Directive 2014/53/EU)

Australian Sponsor:
Emergo Australia, Level 20, Tower II
Darling Park, 201 Sussex Street- Sydney,
NSW 2000
Australia

UK Responsible Person:
Emergo Consulting (UK) Limited c/o Cr360
UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

 MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

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© 2025 Withings. All rights reserved. Apr. 2025